


FDA AUTHORIZED

# Coronavirus (COVID-19) Rapid Test

For Medical Professional Use Only



NOW AVAILABLE: LIMITED SUPPLY

## Our easy-to-use COVID-19 test kit provides results in just 10 minutes



COST-EFFECTIVE



FAST



ACCURATE

### HOW ACCURATE IS THE CORONAVIRUS (COVID-19) RAPID TEST?

Our test has shown in clinical evaluation to have a total agreement of **97.19%** and a kappa value of 0.94.

TEST REAGENT	REFERENCE REAGENT		TOTAL
	POSITIVE	NEGATIVE	
Positive	245	4	249
Negative	16	439	455
Total	261	443	704

Results of the first clinical evaluation, published on March 12, 2020

Clinical evaluation took place across 10 hospitals in China during February 2020 and included 704 tests. Ongoing clinical evaluation and other information is available upon request.

**“I recommend Medek Health Sys. COVID-19 tests 100%, they have a reliable product, excellent customer service, and are very trustworthy and dependable.”**

— Dr. Fernandez, Third Tooth Dental Care, April 16, 2020

### COVID-19 “CORONAVIRUS” IGG/IGM RAPID TEST KIT

To support public health & safety, we are offering a proven, cost-effective option to detect and help stop the spread of coronavirus.

- FDA Authorized (EUA)<sup>1</sup>
- Rapid results within 10 minutes
- Small sample sizes
- Sold in packs of 100
- Following the incubation period, IgM may appear in blood within 3-5 days. IgG will appear as soon as 1-2 weeks.
- Shelf life of up to 24 months from manufacture date
- Forensic / Medical Professional Use Only<sup>2</sup>
- Tests should be conducted by a licensed phlebotomist, or a medical professional
- Test to be performed in a CLIA moderate- or high-complexity setting<sup>3</sup>
- Verification of use case prior to shipping is mandatory

### MORE INFO ABOUT OUR RAPID CORONAVIRUS (COVID-19) TEST KIT

For additional information, such as instructions, how-to videos, clinical trial info, the product insert, and CE documentation, ask your sales representative or visit:

[www.medekhealth.com](http://www.medekhealth.com)

<sup>1</sup>These tests are authorized by the FDA for emergency use as of May 29, 2020. Forensic/medical professional only: Tests should be conducted by a licensed phlebotomist, or a medical professional. Verification of use case prior to shipping is mandatory.

<sup>2</sup>Forensic/medical professional use only: Tests should be conducted by a licensed phlebotomist, or a medical professional. Verification of use case prior to shipping is mandatory.

<sup>3</sup>Please defer to your local regulations. To learn more, download the FDA Fact Sheet: Antibody Test Oversight and Use For COVID-19 on our website.

### AVAILABLE UPON REQUEST:

- Package Insert
- Clinical Evaluation Report
- Stability Report

**“We ran 12 tests on known positives (10) and negatives (2) and got 100pct concordance on the IgG with some of the IgG positives also having varying intensity of IgM. There were no isolated IgM positives but all of these people were sick over 10 days prior and now recovered.”**

— Clinic in East Hampton, NY, April 16, 2020